



## RUDAQUICK<sup>®</sup>-Strips

### Product Description and Purpose

RUDAQUICK<sup>®</sup> plaster strips consist of a transparent, water resistant foil and are provided with a hypoallergenic polyacrylate adhesive. Due to its perforated foil the plaster remains air-permeable. The wound pad is coated so that it does not stick to the wound. Owing to the new individual packing the plaster strip is very easy and quick to apply. By pulling both ends of the packing, one half of the plaster strip is uncovered and can now be fixed at the required spot. The remaining half of the packing is then removed. In this way the plaster strips can be applied very rapidly under hygienic conditions. The strips are used to cover injection sites, small punctures and primary healing smaller injuries.

### Composition

polyethylene vinyl acetate, polyacrylate adhesive

### Normative and Legal Requirements

RUDAQUICK<sup>®</sup>-Strips is a medical product in accordance with the MDD, class I, rule 4, and complies with the requirements of the directive 93/42/EEC.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and EN 980 labels on all its packaging.

### Packaging

Primary packaging: paper packaging

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

### Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

