



## RUDADERM<sup>®</sup>-universal

### Product Description and Purpose

The carrier material is made of a flexible, brown, waterproof synthetic which is coated with a hypoallergenic polyacrylate adhesive. The wound contact layer is protected on all sides by the carrier material, except the smallest format 0,9 x 3,8 cm. The strips are individually packed and are used to cover injection sites, small punctures (for example vein punctions) and primary healing smaller injuries.

### Composition

Polyethylene, polyacrylate adhesive

### Normative and Legal Requirements

RUDADERM<sup>®</sup>-universal is a medical product in accordance with the MDD, class I, rule 4, and complies with the requirements of the directive 93/42/EEC.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and EN 980 labels on all its packaging.

### Packaging

Primary packaging:	paper packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

### Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

