



# NOBA Verbandmittel

Danz GmbH u. Co KG

## NOBAKAWA®

### Product Description and Purpose

The transparent swab is used for cryotherapy after having been cooled in the freezer. The swab's remaining flexible at a temperature of -18°C has proven to be advantageous because it adapts to the body shape. NOBAKAWA® can also be warmed up in the water bath and can be used for thermal therapy. The swab has a protective cover made of white non-woven. NOBAKAWA® can be used repeatedly.

### Composition

Polyethylene, polyamide, water, CMC, propylene glycol, polypropylene

### Normative and Legal Requirements

NOBAKAWA® is a medical product in accordance with the MDD, class I, rule 1, and complies with the requirements of the directive 93/42/EEC.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and EN 980 labels on all its packaging.

### Packaging

Primary packaging: protective cover made of non-woven

Secondary packaging: carton made of cellulose

### Storage

No specific storage conditions to be observed



Date of information: May 18, 2011  
replaced: June 10, 2010