



## RUDAFILM<sup>®</sup>-sterile

### Product Description and Purpose

The carrier fabric of the sterile adhesive dressing consists of a transparent, elastic film made of polyurethane, and is provided with a hypoallergenic polyacrylate adhesive. The wound contact layer, which does not adhere to the wound, is overlapped on all sides by the film. The film protects the wound from the penetration of pathogenic germs. It is waterproof, however steam permeable. The patient is therefore able to manage his daily body care without soaking the dressing. The transparent film allows an assessment of the moisture status of the wound contact layer without removing the dressing. These so-called "island" dressings are individually packed in sterile packaging.

### Composition

Viscose, polyester, polyurethane, polyethylene, polyacrylate adhesive

### Normative and Legal Requirements

RUDAFILM<sup>®</sup>-sterile is a medical product in accordance with the MDD, class Is, rule 4, and complies with the requirements of the directive 93/42/EEC.

Sterilization of the product complies with DIN EN 11137.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and EN 980 labels on all its packaging.

### Packaging

Primary packaging: paper-paper packaging

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

### Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

